

Technical Info



Nutrient Dry Media

Summary acc. to usage for microorganisms

Art. No.	Medium	Media Specialised for...													
		Ph. Eur.	Water	Whole Germ	Coliform	Clo & Bruc & Camp	Enterobac.	Yeast Mould	Lacto & Bac	Pseudo & Leg & Lis	Sal & Shig	Staph & Stre	Chromogen	Special	Mol Biol
X967	2x YT Agar														
X966	2x YT Broth														
6676	2x YT Broth, granulated														
CP61	Acetamide Broth (Solution A)		Water							Pseu					
CL53	Antibiotic Medium No. 1												Special		
CL54	Antibiotic Medium No. 23												Special		
CL79	Azide Dextrose broth						Enterobac.								
X913	Baird Parker Agar											Staph			
AE24	Bengal Red Agar (with chloroamphenicol) (RBC)							Yeast Mould							
CP62	Bile Esculine Azide Agar		Water				Enterobac.								
8281	Blood Agar No. 2 (Base)								Bac	Lis		Strep			
X914	Blood Agar (Base)											Strep		Special	
X915	Brain Heart Infusion Agar			Whole Germ											
X916	Brain Heart Infusion Broth			Whole Germ											
CL06	Brilliant Green Bile Broth		Water		Coliform										
AE43	Brilliant Green Lactose Sucrose Agar										Sal				
AE42	Brilliant Green Tetrathionate Bile Broth										Sal				
5752	Brucella Agar					Bruc									
5754	Brucella Bouillon					Bruc									
5756	Campylobacter Agar, blood-free (Base)					Camp									
8847	Campylobacter Selective Agar (Base)					Camp									
2922	Candida Chromogenic Agar							Yeast Mould					Chromogen		
X937	CASO Agar (Soy Broth, TSB Agar)			Whole Germ											
8848	CASO Agar with Lecithin/Polysorbate 80			Whole Germ										Special	
X938	CASO Bouillon (Soy Broth, TSB)			Whole Germ											
6678	CASO Bouillon, granulated (Soy Broth, TSB)			Whole Germ											
0139	CASO Bouillon, vegetal (Soy Broth, TSB)			Whole Germ											
X918	Cetrimide Agar									Pseu					
2835	CLED Agar													Special	
2836	CLED Agar plus													Special	
2923	Coli Fluoro Medium				Coliform								Chromogen		
CL45	Coliforme Chromogenic Agar		Water		Coliform								Chromogen		
8849	Coliforme Chromogenic Agar (ISO)		Water		Coliform								Chromogen		
X919	Columbia Agar (Base)			Whole Germ										Special	
2997	Cronobacter Chromogenic Agar				Coliform		Enterobac.						Chromogen		



Technical Info

Well advised with Roth.

Art. No.	Medium	Media Specialised for...													
		Ph. Eur.	Water	Whole Germ	Coliform	Clo & Bruc & Camp	Enterobac.	Yeast Mould	Lacto & Bac	Pseudo & Leg & Lis	Sal & Shig	Staph & Stre	Chromogen	Special	Mol Biol
8850	D/E Neutralisierungs-Agar														
AE44	Desoxycholat Citrate Agar										Sal, S hi				
AE25	Dichlorane Bengal Red Agar (with chloroamph.) (DRBC)														
AE26	Dichlorane Glycerol Agar (DG18) (Base)														
CL78	Differential Reinforced Clostridial Medium (DRMC)					Clo									
8295	DNase Test Agar											Staph			
8310	EC Broth														
X920	Endo Agar														
HP17	Fraser Broth (Base)										Lis				
HP18	Fraser Half Broth (Base)										Lis				
HP07	Gelatine Agar (DEV)														
CP64	King's B Medium (Base)										Pseu				
5758	Kligler Iron Agar														
X921	Lactose Broth											Sal			
CL76	Lactose Broth (DEV)														
CL77	Lactose Peptone Broth (DEV)														
AE46	Lactose Sulfite Broth					Clo									
CL80	Laurylsulphat Broth														
X965	LB Agar (Lennox)														
6671	LB Agar (Lennox), granulated														
0159	LB Agar (Lennox), vegetal														
X969	LB Agar (Luria/Miller)														
6675	LB Agar (Luria/Miller), granulated														
X964	LB Broth (Lennox)														
6669	LB Broth (Lennox), granulated														
0155	LB Broth (Lennox), vegetal														
X968	LB Broth (Luria/Miller)														
6673	LB Broth (Luria/Miller), granulated														
CL48	Legionella Agar (Base)										Leg				
5759	Lethen Agar modified														
5761	Lethen Broth modified														
3008	Listeria Chromogenic Agar										Lis				
8322	Lysin Iron Agar										Sal, S hi				
X922	MacConkey Agar														
AE47	MacConkey Broth														
X923	Malt Extract Agar														
CP75	Malt Extract Broth														
CL81	Mannitol Salt Agar											Staph			
5765	Mannitol Egg Yolk Polymyxin Agar (Base)									Bac					
CP73	Marine Broth														
5766	Maximum Recovery Solution														
CL55	m-CP Agar (Base)					Clo									
3014	mEI Chromogenic Agar														
AE65	Mossel-Enrichment Broth										Sal				
X924	MRS Agar									Lac					





Technical Info

Well advised with Roth.

Art. No.	Medium	Media Specialised for...													
		Ph. Eur.	Water	Whole Germ	Coliform	Clo & Bruc & Camp	Enterobac.	Yeast Mould	Lacto & Bac	Pseudo & Leg & Lis	Sal & Shig	Staph & Stre	Chromogen	Special	Mol Biol
HP65	MRS Agar pH 5.7								Lac						
X925	MRS Broth								Lac						
HP64	MRS Broth pH 5.7								Lac						
3026	MRSA Chromogenic Agar											Staph			
1E1K	MR/VP Broth														
X926	Mueller Hinton Agar											Staph			
X927	Mueller Hinton Broth														
8323	Mueller Kauffmann Broth (Base)										Sal				
X928	Nutrient Agar														
CL02	Nutrient Agar (DEV)														
CP66	Nutrient Agar (DEV/ISO)									Pseu					
X929	Nutrient Broth														
X974	NZCYM Broth														
X973	NZYM Broth														
HP09	Orange Serum Agar														
5767	Oxford Listeria Agar									Lis					
5768	Palcam Listeria Agar									Lis					
X917	Peptone Water, buffered														
6680	Peptone Water, buffered, granulated														
AE67	Peptone Water, buffered (Ph. Eur.)														
X930	Plate Count Agar														
HP08	Plate Count Agar with Powdered Milk														
X931	Potato Extract Glucose Agar														
CP74	Potato Extract Glucose Broth														
3075	Pseudomonas Chromogenic Agar									Pseu					
CP67	Pseudomonas Selective Agar (CN Agar) (Base)									Pseu					
CL01	R2A Agar														
CL03	Rappaport Vassiliadis Broth (RVS)										Sal				
AE66	Reinforced Clostridial Medium (RCM)														
8329	Rogosa Agar									Lac					
AE22	Sabouraud 2% Glucose Agar														
AE23	Sabouraud 2% Glucose Broth														
X932	Sabouraud 4% Glucose Agar														
X933	Sabouraud 4% Glucose Agar (with chloramphenicol)														
AE48	Sabouraud 4% Glucose Agar (with chloramph.) (Ph.Eur.)														
CL46	Salmonella Chromogenic Agar										Sal				
X934	Salmonella Shigella Agar										Sal				
5771	Schaedler Agar														
5772	Schaedler Bouillon														
5774	Simmons Citrate Agar														
X935	Slanetz Bartley Agar											Strep			
CP68	Slanetz Bartley Medium (with tetrazolium chloride)														
AE27	SOB Broth														
5778	Sorbitol MacConkey Agar														





Technical Info

Well advised with Roth.

Art. No.	Medium	Media Specialised for...													
		Ph. Eur.	Water	Whole Germ	Coliform	Clo & Bruc & Camp	Enterobac.	Yeast Mould	Lacto & Bac	Pseudo & Leg & Lis	Sal & Shig	Staph & Stre	Chromogen	Special	Mol Biol
AE91	Standard Count Agar														
AE93	Standard Nutrient Agar I														
AE92	Standard Nutrient Broth I														
3935	TAT-Bouillon														
CL47	TBX Chromogenic Agar														
8333	TCBS Agar														
X972	Terrific Broth														
0161	Terrific Broth, vegetale														
HP61	Terrific Broth modified														
3556	Terrific Broth granulated														
1E12	Thioglycolate Broth														
X936	Todd Hewitt Broth														
AE49	Triple Sugar Iron Agar (Base)														
CP69	Tryptone Bile Agar (TBA)														
CP70	Tryptone Soy Agar (TSA)														
CP71	Tryptophane Broth														
CL51	TSC Agar (Base)														
CP65	TTC Chapman Agar (with sod. heptadecylsulph.) (Base)														
8337	Universal-Bier-Agar (Basis)														
3089	UTI Chromogenic Agar														
X939	Violet Red Bile Agar with Dextrose (VRBD)														
AE45	Violet Red Bile Agar with Dextrose and Lactose														
X940	Violet Red Bile Agar with Lactose (VRB)														
CL04	Wort Agar														
5779	Wort Broth														
X941	XLD Agar														
CP63	Yeast Extract Agar														
CL00	Yeast Glucose Agar (with chloramphenicol)														
HP26	Yeast Nitrogen Base (without Amino Acids)														
5780	Yersinia Selective Agar (Base)														
X971	YPD Agar														
X970	YPD Broth														





Well advised with Roth.

Technical Info

Legend:

Ph. Eur.	Media acc. to European Pharmacopeia
Water	Media acc. to DEV or drinking water directives
Whole Germ	Detection of whole germs or unspecific microorganisms
Coliform	Detection of coliforms and <i>E. coli</i>
Clo & Bruc & Camp	Detection of <i>Clostridia</i> , <i>Brucella</i> or <i>Campylobacter</i>
Enterobac.	Detection of enterobacteria
Yeast & Mould	Detection of yeasts or moulds
Lacto & Bac	Detection of lactobacilli or bacilli
Pseudo & Leg & Lis	Detection of pseudomonads, legionells or listeriae
Sal & Shig	Detection of salmonellae or shigellae
Staph & Stre	Detection of staphylococci or streptococci
Chromogen	Chromogenic and fluorogenic media
Special	Media for special applications
Mol Biol	Media for applications in molecular biology

j.h. 09/2020